



Lotus Technologies 2008 Catalogue

Single Use Sterility Test Kits

HYT Single use sterility test kits are used for pharmaceutical product sterility test by the membrane filter method

FEATURES

1. Be manufactured to ISO 9001 and ISO 14001 schedule
2. 28 types of canister are available for different samples
3. New inlet design to eliminate foaming
4. 100% product air testing during manufacturing to ensure product seal integrity
5. Physical tests and Biological test be done before every lot product to be released
6. Perfectly work with HYT Pump and compatible with other sterility test pumps as well



HTY Sterility Test Kits



HTY Sterility Test Kit In Testing

Spcification

Pore size of filter membrane	0.45um HTY membrane
Filter area	14.7 C m ²
Pore size of air filter	0.22um PVDF
Sample container capacity	120ml with graduation marks at 50,75 and 100ml
Tube length	90cm—150ml,depend on types
Max. operating pressure	3 bar (0.3MPa)
Max. operating temperature	60 °C
Sterilization	ETO



HTY Sterility Test Kits with HTY Pump

P O Box 959
Wentworthville
NSW 2145
T: 02 9704 2313
F: 02 9613 0879
sales@lotustech.com.au

www.lotustech.com.au

Lotus Technologies Pty Ltd

